

# PREVENTION (STP)

## TOOTH BRUSHING DEMONSTRATION MODEL



### PE-STP001

#### Features

- Tooth brushing model for individual education and training
- Half-erupted tooth and crowded teeth conditions are represented
- Implant and bridge restorations on the lower right posterior region

#### Purpose

- To demonstrate and instruct brushing techniques for various abnormal teeth alignment conditions
- Teach proper brushing techniques for bridges and implants

## TOOTH BRUSHING DEMONSTRATION MODEL



### PE-STP002 (Permanent dentition)

#### Features

- Adult tooth brushing instruction model with removable teeth (#14 & #16)
- 2 times the natural size
- Lightweight construction
- Dental floss included

#### Purpose

- To demonstrate normal tooth brushing, interdental brushing and dental flossing
- Ideal for group education

\* Tooth numbers indicated above follow the FDI Tooth Numbering System.  
\* Toothbrush sold separately.



Dental floss



Interdental brushing




Brushing



#14 and #16 teeth are removable.

## TOOTH BRUSHING DEMONSTRATION MODEL



### PE-STP004

#### Description

- Toothbrush for PE-STP002, PE-STP003

### PE-STP005

#### Description

- Interdental brush for PE-STP002, PE-STP003

## TOOTH BRUSHING DEMONSTRATION MODEL



### PE-STP003

#### Features

- Pediatric tooth brushing instruction model with primary dentition
- 2.5 times the natural size
- Lightweight construction

#### Purpose

- To demonstrate normal tooth brushing, interdental brushing and dental flossing
- Ideal for group education

\* Toothbrush sold separately.

## CARIES PREVENTION MODEL



### PE-STP007

#### Features

- 2x size model
- Upper posterior teeth show caries metal inlay restoration as well as sealant for caries prevention

#### Purpose

- For explaining sealant application to developmental groove areas for the purpose of caries prevention
- For explaining progressed caries and treatment by metal inlay restoration

## CARIES PREVENTION MODEL



### PE-STP008

#### Features

- Various caries conditions shown on the upper right premolar and molar
- Sealant treatment demonstrated on the healthy second molar

#### Purpose

- Explain the use and effectiveness of sealants for caries treatment and prevention

## PROSTHESIS CARE AND BRUSHING MODEL



### PE-STP011

#### Features

- Removable partial denture and 3-unit implant bridge
- 2 times the natural size with tongue
- Lightweight construction

#### Purpose

- Proper brushing demonstration to patients that have partial dentures and implant restorations